

Date: September 8, 2008

Logger: JL,AB

Sample Number	From (m)	To (m)	Interval (m)	Assay REC		Weight (kg)	Note: Standard Number, Blanks, Dup's
				Metres	Percent		
ML03478	2.30	5.00	2.70	2.70	100	5	
ML03479	5.00	8.00	3.00	3.00	100	4.5	
ML03480	8.00	11.00	3.00	3.00	100	4.75	
ML03481	11.00	13.95	2.95	2.92	99	4.5	
ML03482	13.95	17.00	3.05	2.97	97	5.5	
ML03483	17.00	20.14	3.14	3.12	99	5.75	
ML03536	20.14	23.00	2.86	2.86	100	5.75	
ML03537						3	BLANK
ML03538	23.00	26.00	3.00	2.94	98	5.5	
ML03539	26.00	29.00	3.00	3.00	100	5.5	
ML03540	29.00	31.98	2.98	2.88	97	5.25	
ML03541	31.98	35.00	3.02	3.00	99	5	
ML03542	35.00	38.31	3.31	3.26	98	6	
ML03543	35.00	38.31	3.31	3.26	98	2.3	DUPLICATE
ML03544	38.31	41.00	2.69	2.69	100	4.5	
ML03545	41.00	44.00	3.00	3.00	100	5	
ML03546	44.00	46.50	2.50	2.50	100	4.75	
ML03547							STANDARD ML-5
ML03548	46.50	48.95	2.45	2.45	100	4.5	
ML03549	48.95	51.22	2.27	2.27	100	3.75	
ML03550	51.22	53.30	2.08	2.08	100	4	
ML03551	53.30	55.50	2.20	2.17	99	3.5	
ML03552	55.50	58.00	2.50	2.45	98	4	
ML03553	58.00	61.12	3.12	3.09	99	4.4	
ML03554	61.12	64.00	2.88	2.81	98	5.5	
ML03555	64.00	66.78	2.78	2.78	100	5.5	
ML03556	66.78	69.49	2.71	2.67	99	4.2	
ML03557	69.49	71.91	2.42	2.42	100	4.1	
ML03558						3	BLANK
ML03559	71.91	75.00	3.09	3.05	99	5.4	
ML03596	128.75	131.47	2.72	2.72	100	4.5	
ML03597	128.75	131.47	2.72	2.72	100	2	DUPLICATE
ML03598	131.47	134.50	3.03	2.93	97	5	
ML03599	134.50	137.50	3.00	3.00	100	5	
ML03600	137.50	139.94	2.44	2.25	92	4	
ML03601	139.94	142.84	2.90	2.88	99	6.1	
ML03602	142.84	145.23	2.39	2.35	98	4.75	
ML03603						3	BLANK
ML03604	145.23	148.25	3.02	3.02	100	6.4	
ML03605	148.25	151.31	3.06	3.01	98	5.5	
ML03606							STANDARD ML-5
ML03607	151.31	154.00	2.69	2.69	100	5.5	
ML03608	154.00	156.34	2.34	2.34	100	2.8	
ML03609	156.34	158.00	1.66	1.58	95	3.25	
ML03610	158.00	160.08	2.08	2.05	99	4.3	
ML03611	160.08	162.56	2.48	2.47	100	4	
E.O.H.							